

# High resolution ultrasound scanner for skin imaging

Christine Turlat Sales Director Atys medical 17 Parc d'Arbora 69510 SOUCIEU EN JARREST



High resolution ultrasound scanner

Atys company Principle of ultrasound imaging DERMCUP Normal image of skin Examples of applications Conclusion



# Atys medical



17 Parc d'Arbora 69510 SOUCIEU EN JARREST NEAR LYON FRANCE Year of creation: 1990

Privately owned company

Activity: design, manufacturing & sales of medical devices

French leader in its field



High resolution ultrasound scanner

# ATYS medical range of products





# Principle of ultrasound

- Medical ultrasound imaging is ultrasound that is converted to an image
- Diagnostic Medical applications in use since late 1950's
- Frequency ranges used in medical ultrasound imaging are 2 - 15 MHz General abdomen, OB/Gyn: 3.5MHz
- The use of these frequencies for the diagnosis of process (size of 0.1 mm) limited to the corium is inadequate or impossible. They do not provide satisfactory image of the dermis.
- DERMCUP: 25 MHz





High resolution ultrasound scanner

# Principle of ultrasound

### **Transducer frequency**

#### Higher frequency = higher axial resolution Lateral resolution depends on the beam Higher Frequency = lower penetration

Resolution: ability to distinguish two structures that are close together as separate.

# The ideal situation is to use the highest frequency possible to achieve penetration to the area of interest.

Typical diagnostic ultrasound for fetal imaging is around 3.5 to 7 MHz. DERMCUP uses 25 MHz central frequency probes. This makes it suitable for imaging the skin and superficial soft tissue to a depth of approximately 1 cm.







# Principle of ultrasound

The probe includes a piezoelectric crystal The sound waves are focused into a beam and transmitted in the soft tissues of the body.

The transducer pauses to received the reflected waves

- When the sound waves hit a boundary between acoustically different tissues some of the waves are reflected back.
- Other waves travel further until they reach another boundary.
- The reflected waves interact with the piezoelectric crystal to produce a small current. These signals are processed by the ultrasound machine that calculates the location and intensity of these reflections.







High resolution ultrasound scanner

#### How is the B-mode format image formed on the monitor?

- •The imager displays a real time two dimensional image.
- The strength or amplitude (brightness) of each reflected wave is represented by a dot.
  The position of the dot represents the depth from which the returning echo was received
  These dots are combined to form a complete image
- •The DERMCUP uses a single element transducer driven by a small motor that transverses back and forth within the probe.







# **DERMCUP** features

Frequency:16, 20, 25 MHz

### Sectorial probe

 Scanning length: 6 mm produces sector or pie-shaped image

#### Linear probe

- Scanning length: 16 mm

#### Resolution @ 25 MHz

- axial: 60 μm
- lateral: 150 μm

Exploration depth: 12 mm







## High resolution ultrasound scanner

### DERMCUP probe Easy to use

Probe perpendicular to the skinApply minimum pressureApply gel between probe and skin

# Easy to maintain

-Water chamber filling -Membrane replacement











Horizontal structure of the skin or fibrillar network comprised of collagen and elastic fibers. This is the fribrillar network that is responsible for ultrasound echogenicity of the dermis.

The dermis provides a natural contrast medium in which different pathologies can be outlined if they cause low reflectancy or disturbance of interface and dimensions.

Atys<br/>medicalHigh resolution ultrasound scannerDERMCUP skin image

Structures such as blood vessels, glands, fascia and adipose tissue can be seen in the hypodermis.
The deep fascia can be identified below a layer of subcutaneous tissue
Muscles have few internal echoes.
Bones cause heavy reflection.









High resolution ultrasound scanner

The thickness of dermis varies with anatomic area.





# **DERMCUP** applications

- · Vertical skin cross section in vivo
- · Accurate non invasive measuring device for skin pathology
- Skin cancer management applications
  - Assessment of the size of the tumours to aid in surgical planning
  - Improvement of the accuracy of the clinical diagnosis
  - Assessment of the efficiency of the treatment
- · Objective 2 D documentation for most treatment modalities
- · Quantitative 2 D evaluation of response to topical treatments





Pre-op Breslow's 3.68mm

(histology 3.7mm Level IV)

Recurrence nBCC on left nostril (1 year after surgery)





## High resolution ultrasound scanner



#### Histiocytofibroma



#### Seborrhoeic keratosis



#### **Epidermal cyst**





Subepidermal bubble





Angioma of the lip

Seborrhoeic wart





psoriasis



Helpful in diagnosis of glomus tumors which are often located under the nail





Normal finger tendon (longitudinal)





Tendon with inflammation (longitudinal & transversal)



## High resolution ultrasound scanner

Gum





# Eye



Eye



Injured eye





### **Dermo cosmetic applications**

•Quantitative and objective evaluation of treatments

•Thickness of the dermis (Ultrasound imaging demonstration of the improvement of non-ablative laser remodeling by concomitant daily topical application of 0.05% retinaldehyde)

•Collagen concentration: brightness

•Ultrasound reveals the appearance of a subepidermal low echogenic band that thickens with age, especially in environmentally exposed areas.



Atys medical

High resolution ultrasound scanner

### Advantages of ultrasound imaging

- Non invasive
- · Safe for the patients and the users
- Can be repeated as often as necessary
- Quite affordable
- Quantitative and objective evaluation
- Not time consuming
- No special preparation of the patient



# **FUTURE EVOLUTIONS**

- Higher frequency: up to 50 MHz
- Implementation of a Doppler module to study the vascular flow
- 3 D probe





High resolution ultrasound scanner

# Thank you for your attention