







Introduction

New techniques...why?

- √ Interest of new technologies for the imaging of the skin
- ✓ New & mature Technologies "Research and routine"
- ✓ Imaging "to perceive" & imaging "to understand"
- √ Steps of validation of techniques and related tools (contention systems and software)
- √ Tools for standardization and software for quantification









Interest of new technologies for the imaging of the skin

- For the market :
 - √ Needs of innovation in dermocosmetology to support innovation of products
 - √ Needs of « images » to support « science » and to « show »
- For the science:
 - √ More and more complex activity of the skin care product
 - √ Real biological activity



Increase interest of new technologies













New and mature techniques Research vs Routine

- New techiques for research:
 - √ In a specific context
 - √ Involve in a specific environment
 - √ Need large competencies for development of techniques
- Mature techniques for routine:
 - ✓ Reliable, stable, accurate...
 - √ Validated
 - ✓ Related tools for quantification in stable conditions
 - ✓ For "intensive and reliable production of data"



A gap between 2 fields of application
To different use











Imaging « to perceive » Imaging « to understand »

- To perceive:
 - ✓ Quantitative digital photography
 - ✓ Dermoscopy
 - √ Realistic 3D acquisitions
 - √ ...
- To understand
 - ✓ Intra-cutaneous imaging
 - √ Ultrasound scanning
 - ✓ OCT
 - √ Confocal microscopy
 - √ Multiphotonic mocroscopy
 - √ RAMAN spectrocopy
 - V ...

- Objective validation and measurement of the perception (in parallel of the clinical evaluation)
- Standardization & calibration



- To know the mechanism of the product action
- To prove the real activity
 To quantify the target of its activity (and not a consequence)















Imaging « to perceive » Imaging « to understand »

- To perceive:
 - ✓ Quantitative digital photography
 - ✓ Dermoscopy
 - ✓ Realistic 3D acquisitions
 - ✓ ...
- To understand
 - ✓ Intra-cutaneous imaging
 - ✓ Ultrasound scanning
 - **√** OCT
 - √ Confocal microscopy
 - ✓ Multiphotonic mocroscopy
 - √ RAMAN spectrocopy
 - √ ...



- Objective validation and measurement of the perception (in parallel of the clinical evaluation)
- Standardization & calibration




- To know the mechanism of the product action To prove the real activity
- To quantify the target of its
- activity (and not a consequence)

















New technologies transfer ...our job!

Consultancy, Expertise, Training on New Techniques/Methods and Special Protocols

ORION Concept

Laboratory for the development, the validation and use of new methods and techniques for the study of the skin and the efficacy of skin care product (High technologies for imaging)

Collaborations and partnerships
Validation of techniques & Methods
Validation of products with increase value techiques
Skin research



COSMI'O

Sales of « Packages » (Matrials, sofwares, Training&Support)

- · Materials for data acquisitions
- Validated sofwares for quantifications
- · Installations, Training and Support for users...













Exemple 1 : Quantitative digital photography of the skin : easy?

- Definition of needs:
 - √ Color (Product or skin)
 - √ Homogeneity (Skin complexion, make-up products, ageing....)
 - ✓ Detection and quantification of different elements onto the skin (Pigmentary spots, telangiectasies, wrinkles and their density, acne lesions, hairs ...)
 - Measurements of lashes (Real and apparent length, width, volume and curvature)
 - √ Visibility of wrinkles (On pseudo photographs)
 - √ Etc...


On parametric recontructed images (Spectroscopic images)


 Selective quantifications on only one component of the skin (Vascularisation or pigmentation)


















Software for the controlled quantifications


FrameScan_®V03








Red Channel



RGB channel separation *

Green Channel



Blue Channel

Red: High homogeneity of the skin - Usefull for lashes and hair evaluations

Bleu: High contrast for pigmentary and vascular elements onto the skin















Software for the controlled quantifications

Reconstructed parametric images





Acquired Image (Crossed polarized image)



Pigmentation



Vascularisation



Contrast

Separation of different components of the skin


QUANTIFICATIONS:


- Intensity, size, homogeneity, distribution, density (Acne, telangiectasies, complexion, pigmentation,
- Multi-area measurements (Local & Global cartography)
- "Contrast of the face" Perception of the relief and ptosis (Cutaneous slackening)

















Software for the controlled quantifications

FrameScan_®V03











Acquired image (crossed polarized)

Red channel

Perceived Lenght Without product

Real length Without product

Real length WITH product

QUANTIFICATIONS: With "active contours"

- Real & Apparent Length, width (Validated correlation with the clinical evaluation)
 -Number, crossings (parallelism of lashes) , aggregates (resolution <8µm)

And a lot of other functions...













Exemple 2 : Quantification of Confocal Laser Microscopy ConfoScane V01.1.0

A mature technique and system for routine

A necessary extraction of quantitative information from the images

A large amount of data for analysis requiring a specific management



Management of data and sample images (in stacks or blocks)

Multi-Selection ROIs for targeted analysis

ORION Concept





