



**IN VIVO
SKIN STUDIES BY
RAMAN
SPECTROSCOPY**

Johanna de Sterke, André van der Pol
 Peter Caspers, Rolf Wolthuis
 River Diagnostics B.V.
 The Netherlands



River Diagnostics - *in vivo* skin analysis

- River Diagnostics founded in 2002 as a spin off from the Erasmus University Medical Centre Rotterdam
- Model 3510 Skin Composition Analyzer in the commercial market since 2004



C.E.R.I.E.S.
 ESTÉE LAUDER

Kimberly-Clark

SHISEIDO

proDERM

Johnson & Johnson

COLGATE-PALMOLIVE COMPANY

Kanebo

SSL International plc

P&G






Principles of measurement - Raman skin analysis

- optical, non-invasive and non-disturbing
- no biopsies
- objective
- high spatial resolution
- no dyes or labels involved
- detailed information about chemical composition




Raman spectroscopy - principle





Information in Raman spectra of skin

5




It is all about:



Main chemical composition of the skin at different depths



confocal Raman microspectroscopy

depth resolution: ~5 µm







Information in Raman spectra of skin, examples

9

- skin characterisation
 - water, nmf, lipids, carotenoids, urea
 - typing, aging effects, ethnic effect, cultural effects
 - neutraceutical
- Penetration and penetration kinetics
 - glycerol, dmso, lactate, ethanol, caffeine, menthol, nicotinamide, uv absorbers, oils, retinol, carotenes
 - formulation optimization
 - transdermal drug delivery, various pain relief agents
 - Toxicological: pyrene, toluene, butoxyethanol
- metabolic processes
 - hydrolysis of esters
- effects of products
 - moisturization
 - pH modulation
- effects of environment / skin treatment
 - humidity, season, climate, solar exposure
 - bathing, washing
- medical (diseased/disfunctional skin)
 - Atopic dermatitis



Applications

10

- Water concentration measurement
- Skin characterization
- Penetration of active






applications – effects of moisturizers

In 2006, a large study was published by PG researchers:

Sieg A, Crowther J, Blenkiron P, Marcott C, Matts PJ. Confocal Raman Microspectroscopy – Measuring the effects of topical moisturizers on stratum corneum water gradient in vivo. Proc. of SPIE Vol.6093 60930N, 2006.





applications – effects of moisturizers





applications – effects of moisturizers





River
Diagnostics



skin hydration: swelling

14



River
Diagnostics



Applications


15

- Water concentration measurement
- Skin characterization
- Penetration of active



applications – characterization of aging

urea





Applications

17

- Water concentration measurement
- Skin characterization
- Penetration of active





River
Diagnostics



How to extract practical information, concentration profiles

18




River
Diagnostics




How to extract practical information, kinetics of diffusion

19




How to extract practical information, Absorbed volume of product

20



Absorbed volume product



Ingredient Concentration in mmol/gr keratin \rightarrow Absorbed amount in $\mu\text{g}/\text{cm}^2$





applications - retinol delivery

21

Pudney et al., Applied Spectroscopy 61(8): 2007

An *In Vivo* Confocal Raman Study of the Delivery of *Trans*-Retinol to the Skin

PAUL D. A. PUDNEY,* MICKAËL MÉLOT, PETER J. CASPERS, ANDRE VAN DER POL, and GERWIN J. PUPPELS

Measurement Science Unit, Unilever R&D, Colworth Science Park, Sharnbrook, Bedford, MK44 1LQ (P.D.A.P., M.M.); and River Diagnostics B.V., Rotterdam, The Netherlands (P.J.C., A.v.D.P., G.J.P.)

- **Solution 1:** 70 % ethanol, 30 % PG and 0.3 % *trans*-retinol
- **Solution 2:** 99.7 % MYRITOL®318 and 0.3 % *trans*-retinol
- monitor penetration for 0-6 hrs after application




applications - retinol delivery

22

Pudney et al., Applied Spectroscopy 61(8): 2007

in PG/EtOH



in MYRITOL®318

